



**CPAL**

Central Pennsylvania Alliance  
Laboratory

# Technical Bulletin

**No. 166**

**April 30, 2018**

## **Hepatitis B Surface Antibody -Assay Restandardization-**

**Contact:** Steph Frey, 717.812.4768  
*Operations Manager, Clinical Pathology, CPAL*

Dr. Jennifer Thebo, 717.851.3210  
*Director, Technical Operations and Scientific Affairs, CPAL*

### **Ordering Information and Suggested Codes:**

<b>Mnemonic:</b>	<b>HBSAB QNT</b>	<b>.HEPPOSTVAC</b>	<b>EHBSAB</b>
<b>Test Name:</b>	Hepatitis B: HBsAb Quant (Total)	Hepatitis B: Post Hepatitis B Vaccination (LGH)	Hbsab Quant (Emp LGH)
<b>Test Number:</b>	1750042	1750043	1750039
<b>Specimen:</b>	Serum	Serum	Serum
<b>LOINC Codes</b>	5193-8	5193-8	5193-8
<b>CPT Codes</b>	86706	86706	86706

**Effective Date:** April 30, 2018

**Performed:** Monday through Saturday

### **Introduction:**

Abbott Diagnostics has enhanced their Hepatitis B Surface Antibody (AUSAB) assay as part of their ongoing commitment to product improvements. The Architect AUSAB assay has been re-standardized to the WHO 2<sup>nd</sup> International Reference Preparation, 2008 (Code 07/164) instead of the WHO 1<sup>st</sup> International Reference Preparation, 1977 (Code W1042). An upward shift of approximately 6% - 14% may be observed with patient specimens upon implementation of the newly re-standardized assay. Abbott conducted studies on negative and positive patient populations. The studies demonstrated that the performance of the 1<sup>st</sup> WHO standardized and 2<sup>nd</sup> WHO standardized assay is in agreement.

### **Validation Studies:**

#### **Precision:**

Two levels of control were run ten times each, within the same run and on two different days. The manufacturer's claim for within run precision of  $\leq 10\%$  for specimens with a positive result was met. Abbott does not publish criteria for specimens with a result of negative. (Table 1)


**Table 1: Precision**

Test	Within Run Precision				Between Run Precision				Accept
	Neg Pool Mean	%CV	Pos Pool Mean	%CV	Neg Pool Mean	%CV	Pos Pool Mean	%CV	
AUSAB	2.70	10.46	896.56	2.20	2.70	10.35	892.73	1.84	Yes

**AMR and Linearity:**

To verify the AMR and linearity of the assay, AUSAB calibrators of a different lot than those used to calibrate the assay, were run in triplicate. The results were entered into EP Evaluator. Please see the EP Evaluator results below (Figure 1). The calculated slope is 0.991 with an intercept of 0.030 and an observed error of 0.043 mIU/mL or 0.9 %. The assay is linear within Allowable Systematic Error of 1 mIU/mL or 20.0%. The linearity validation is acceptable and the AMR of the assay verified.

**Figure 1: AMR and Linearity**



**Method Comparison:**

Forty-three specimens were split and processed utilizing the current Architect AUSAB and the 2<sup>nd</sup> WHO standardized Architect AUSAB assay. Quantitative analysis yielded a correlation coefficient (R) of 0.9955 with a slope of 1.118 and an intercept of 3.033 (Figure 2). Overall bias for the current assay vs the newly standardized assay was 14.997%.

Figure 2: Method Comparison

# EP Evaluator®


Clinical Laboratory – Wellspan Health

AUSAB

## Alternate (Quantitative) Method Comparison

X Method AUSAB

Y Method AUSAB Di



### Regression Analysis

	Deming	Regular
Slope	1.118 (1.084 to 1.151)	1.112 (1.079 to 1.148)
Intercept	3.033 (-4.375 to 10.441)	3.562 (-3.835 to 10.960)
Std Err Est	21.742	21.712

95% Confidence Intervals are shown in parentheses

### Supporting Statistics

Corr Coef (R)	0.9955	SubRange Bounds	None
Bias	14.199 (14.997 %)	Points (Plotted/Total)	43/43
X Mean ± SD	94.676 ± 202.500	Outliers	Not Tested
Y Mean ± SD	108.874 ± 226.269	Scatter Plot Bounds	None
Std Dev Diff	31.269		

### References:

1. Product Information Letter, Abbott Diagnostics, January 2018.
2. Architect AUSAB IFU, Abbott Diagnostics, March 2017.